

Second Harmonic Generation (SHG)/Two-photon Excitation Fluorescence (TPEF) Microscopy Imaging and Quantitative Assessment of Septal Fibrosis in Choline Deficient High Fat Diet (CDHFD) Rat Model of Nonalcoholic Steatohepatitis (NASH)

Dr. Xiao Teng^{†1}, Dr. Rohan Manohar^{†2}, Dr. Hiroaki Yashiro², Qiang Yang³, Alvin Leong¹ and Dr. Gideon Ho¹

(1) Histoindex Pte Ltd, Singapore (HistoIndex), (2) Takeda Pharmaceuticals, Gastroenterology Drug Discovery Unit, Cambridge, USA, (3) Hangzhou Choutu Technology Co., Ltd.

† Both authors contributed equally to this work.



HISTOINDEX®
New Standard | New Life

Animal models of NASH are crucial in elucidating NASH pathogenesis and examining the therapeutic effects of various agents. However, there are very few preclinical models that exhibit features of late-stage fibrosis within a short time, including septae formation that is typical in F3 noncirrhotic and F4 cirrhotic human liver. In addition, there are few quantitative methods to assess the changes of septal fibrosis.

The purpose of this study was to develop an automated quantitative system for septa fibrosis and evaluate it in a CDHFD preclinical rat model of NASH that mimics histological changes in F3 noncirrhotic patients.

EXPERIMENT OVERVIEW

Wistar Han rats were fed CDHFD (Research Diet#A16092003), or control diet (Research Diet#A12450K) for 12 weeks as previously described⁴. The Acc inhibitor (ND630) was dosed from week 6 onwards in 2 separate studies. Liver tissue specimens were processed into formalin fixed paraffin embedded (FFPE) slides.

Stain-free SHG/TPEF microscopy (Genesis®200, HistoIndex, Singapore) was used for imaging, while an artificial intelligence (AI)-based system analyzed septa collagen in the unstained tissues. The algorithm first recognized the portal tract (PT) and central vein (CV) zones, and then examined each bundle of collagen in between using histological diagnosis of septa in human as a gold standard reference. Thereafter, the septae formation and reduction in these studies were detected, quantified and compared among different groups.

FFPE slides were stained with Picrosirius Red (PSR), whole slides were scanned using Leica Versa 200. PSR+ percentage analysis was performed in Leica Aperio ImageScope using the Color Deconvolution module.

Figure 1: Correlation between PSR+ and %SHG area





Figure 2: Septal area identification



1. In both studies, we observed a significant increase in the %SHG measured collagen area and %PSR in the CDHFD Vehicle rats compared to control diets.
2. CDHFD vehicle animals in both studies exhibited significant of septae between PT and CV zones, though the majority were delicate ones that are usually observed as regressive septae in human F3 patients. Control diet fed animals did not exhibit septae formation.
3. ND630 significantly reduced PSR+ area and septal collagen area and showed a trend reduction in %SHG in both studies.

Septae formation is observed in the CDHFD rat model after 12 weeks on diet, and this model might be suitable for septal assessment of NASH and changes with therapeutic agents. Septae in CDHFD rats at the 12th week on diet resemble the regressive septae in human F3 patients. Longer diet induction may be necessary to observe progressive septae. Further evaluation will be necessary to elucidate differential effects of ND630 on %SHG and %Septae.

An automated AI-based stain-free imaging system with SHG/TPEF could be useful for monitoring liver fibrosis activity in a fully quantitative manner, especially septal changes in animal models of NASH.

1) Wanxin Sun et al. "Nonlinear optical microscopy: use of second harmonic generation and two-photon microscopy for automated quantitative liver fibrosis studies." Journal of Biomedical Optics 2008.
2) Feng Liu, et al. "Automated evaluation of liver fibrosis in thioacetamide, carbon tetrachloride, and bile duct ligation rodent models using second-harmonic generation/two-photon excited fluorescence microscopy." Laboratory Investigation 2016.
3) Feng Liu, et al. "gFIBS: A Novel Automated Technique for Quantitative Evaluation of Fibrosis, Inflammation, Ballooning, and Steatosis in Patients with Non-Alcoholic Steatohepatitis." Hepatology 2020.
4) Bates et al. "Acetyl-CoA carboxylase inhibition disrupts metabolic reprogramming during hepatic stellate cell activation." J Hepatology 2020.

This study was in part funded as a fee for service study by Takeda Pharmaceuticals, paid to HistoIndex.

R.M and H.Y are employees of Takeda Pharmaceuticals and own stock in Takeda. X.T, Q.Y, A.L and G.H are employees of Histoindex or its subsidiary, X.T holds stock options and G.H owns stock in HistoIndex.

Dr. Rohan Manohar, rohan.manohar@takeda.com
Mr. Alvin Leong, alvin.leong@histoindex.com